## ENGLISH

### ELECTRONIC PERSONAL SCALE

#### EXCESS WEIGHT SCALES ANALYZER

This appliance is to help you turning to healthy life-style. Read this manual carefully and keep it for future reference.  
The «BIA» technique (Bioelectrical Impedance) (complete body resistance) Analysis implemented by the appliance allows calculating fat, water and muscle mass percentage with high accuracy considering personal data (height, age, gender) stored in memory.

#### WHAT IS EXCESS WEIGHT?

It is well known that excess weight may affect health adversely. Excess weight may cause fatness, heart problems, progress of atherosclerosis, osteoporosis, although the lack of weight may lead to osteoporosis too. (Osteoporosis is the systemic skeleton disease. Its features are as follows: decrease of bone mass per volume, and impairment of bone tissue micro architecture which lead to bones fragility, risk of fractures because of minimal trauma or even without one.)

The physical exercises and moderate diet help to get rid of excess weight, and excess weight analyzer is the ideal tool to control excess weight.  
It is well known, that in case of insufficient weight muscles and tissues contain more water than fat. It is also common that water conducts electric current better than fat does. You can easily measure your body fat percentage in muscles by conducting weak electric current through the body. This helps you keeping an eye on your forms.

#### Recommendations

- Weigh yourself at the same time of a day.
- Remember that 2-3 hours should pass after the last food intake.

Attention! Data on excess weight could be obtained only if you stand on the scales barefoot.

#### Safety Measures

- The pregnant women are ABSOLUTELY prohibited from using the excess weight analyzer.
- The persons having got pacemaker and other medical equipment with built-in accumulators are ABSOLUTELY prohibited from using the excess weight analyzer.

#### Attention!

- In case of misuse of the appliance or its overload, the symbol «Err» (error) appears on the display.
- Symbol «Lo» reminds that the batteries should be replaced.
- While measuring percentage of fat ("FAT"), water ("WATER") and muscle mass ("MUSCLE MASS") symbols «- - -» identify that there is a lack of contact with sensors. It is required to stand on the scales barefoot thus your feet and heels are on the metal plates.

#### Using the scales

Open batteries compartment cover and make sure that 2 lithium CR2-032 batteries are installed properly. Remove the insulating gap to connect the batteries.

#### Selecting units

There is a switch on the lower part of the scales' housing. When it is in "Kg" position, weight would be measured in kilograms and the high would be measured in centimeters. When the switch is in "Lb" position, weight would be measured in pounds and the high would be measured in feet and inches. When the switch is in "St" position, all weight measurement results would be displayed in stones.

#### Switching on

Slightly knock the housing to switch the scales on.

#### FUNCTION BUTTONS

- 1. - decreasing displayed number
- 2. **SET** - selection and confirmation of data
- 3. - increasing displayed number

#### MEASURING WEIGHT

1. Switch the appliance by a slight knock on its housing.
2. First all icons are displayed, then symbols "0.0 KG", "0.0 LB" or "0.0 ST" are displayed depending on units selected by the switch located on the lower part of scale's housing.
3. Stand on the scales. Stay motionless and wait until your exact weight appears on the screen. Moving up symbols «- -» are displayed while weight measuring is going on. Your weight is displayed for a few seconds; afterwards the scales are switched off automatically.

#### Inputting personal data

1. Press button **SET**, user number will be displayed. Increase the number by pressing button «-» or decrease it by pressing button «+». You can input personal data for 12 users

2. Press button **SET** to confirm selected user number.
3. Afterwards, the gender indicator is lighted up:
  - ♂ - for man and ♀ - for woman.
  - ♂ - for men doing sports regularly (10 hours a week).
  - ♀ - for women doing sports regularly (10 hours a week). Press button «+» or «-» to change current symbol. Confirm selection by pressing button **SET**.
4. The appliance is ready for your age data input. Select your age by pressing «+» or «-» buttons. The appliance supports age range from 10 to 99 years old. Confirm selection by pressing button **SET**.
5. The appliance is ready for your high data input. Select your high by pressing «+» or «-» buttons. The appliance supports high range from 100 to 220 cm (11-22 m) or from 3' 03" to 7'02" (-' feet"/ -" inch). Confirm selection by pressing button **SET**.
6. Your personal data are displayed, afterwards the appliance switches to weight measuring mode.

7. Stand on the scales, stay motionless and wait until an indication of your exact weight appears on the display. Moving up symbols «- - -» are displayed while measurement analyzer is going on. Your weight is displayed for several seconds; afterwards the appliance is switched off automatically.

#### Note

1. If time range between each personal data input is more then 10 sec, the appliance is switched to weight measuring mode automatically.
2. A user should to undertake all the actions indicated above (1-4) in order to get the percentage of fat, water and muscle mass.
3. If no any information was provided, the appliance would measure percentage of fat, water and muscle mass for 30 years old man 1,76 cm high.

#### EXCESS WEIGHT MEASUREMENT

1. Select the user (01, 02, ...12) by pressing «+» or «-» button.
2. After selection of user, symbol "0.0 Kg" and personal data (gender, age, high) are displayed on the display.
3. Stand on the scales barefoot, your foot-plates and heels should be fitted into the metal plates. Wait until your weight is displayed.
4. Keep staying on the scales. After some time "0.00" symbols are displayed and begin disappearing gradually (as percentage of fat, water of body tissues is being measured).
5. After several seconds your weight is displayed, afterwards the percentage of fat ("FAT"), water ("WATER") and muscle mass ("MUSCLE MASS") appears.
6. After 16 seconds the appliance will be switched off automatically.

#### Notice:

- In case of unstable performance of the appliance, examine the condition of batteries. If the batteries are discharged, replace them.
- Take care using the scales. The scales are a precision tool; avoid temperature differences, excess humidity, strikes, dust and exposition to direct sun rays.
- Keep away from heaters.
- Do not jump on the scales.
- Wipe the platform with damp cloth.
- Clean the basement of the scales with soft cloth.
- Do not use chemicals and abrasive agents to clean the scales.
- Do not let liquids inside the housing of the scales.
- Do not use or keep the scales inside excessively humid (more than 80%) cabinets.
- Do not allow the appliance getting into a contact with water or other liquids, as such contacts could affect the measurement accuracy or lead to eternal damage of the scales.

#### MAINTENANCE

- Take care using the scales. The scales are a precision tool; avoid temperature differences, excess humidity, strikes, dust and exposition to direct sun rays.
- Keep away from heaters.
- Do not jump on the scales.
- Wipe the platform with damp cloth.
- Clean the basement of the scales with soft cloth.
- Do not use chemicals and abrasive agents to clean the scales.
- Do not let liquids inside the housing of the scales.
- Do not use or keep the scales inside excessively humid (more than 80%) cabinets.
- Do not allow the appliance getting into a contact with water or other liquids, as such contacts could affect the measurement accuracy or lead to eternal damage of the scales.

#### SWITCHING ON

Slightly knock the housing to switch the scales on.

#### FUNCTION BUTTONS

- 1. - decreasing displayed number
- 2. **SET** - selection and confirmation of data
- 3. - increasing displayed number

#### MEASURING WEIGHT

1. Switch the appliance by a slight knock on its housing.
2. First all icons are displayed, then symbols "0.0 KG", "0.0 LB" or "0.0 ST" are displayed depending on units selected by the switch located on the lower part of scale's housing.
3. Stand on the scales. Stay motionless and wait until your exact weight appears on the screen. Moving up symbols «- -» are displayed while weight measuring is going on. Your weight is displayed for a few seconds; afterwards the scales are switched off automatically.

This product conforms to the EMC-Requirements as laid down by the Council Directive 89/336/EEC and to the Low Voltage Regulation (73/23 EEC).

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